**Introductie**

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BC-37, Non-contact Infrared Forehead / Object Thermometer

**System Owner**

Thank you for purchasing the Universal Infrared Forehead Thermometer. Please read this instruction manual first so you can use this thermometer safely and correctly. Please keep this instruction manual for future reference. This innovative medical device uses advanced infrared (IR) technology to measure temperature instantly and accurately on the forehead or object.

Always check our website for the latest version of this product's user instructions.

**IMPORTANT SAFETY INSTRUCTIONS, READ BEFORE USE**

- The following basic safety precautions should always be taken:
1. Close supervision is necessary when the thermometer is used by, on, or near children, handicapped persons or invalids.
2. Use the thermometer only for the intended use described in this manual.
3. Do not use the thermometer if it is not working properly, or if it has suffered any damage.

**BEFORE YOU BEGIN**

**Cautions and Warnings:**

- As with any thermometer, proper technique is crucial to obtaining accurate temperature readings. Please read this manual thoroughly and carefully before using.
• Always operate the thermometer in an operating temperature range 16°C to 40°C (60.8°F to 104°F), and relative humidity less than 85%.
• Always store the thermometer in a cool and dry place -25°C to 55°C (-13°F to 131°F) and relative humidity less than 95%. Avoid direct sunlight.
• Avoid dropping the thermometer, if it happens and you think the thermometer may be damaged, please contact customer services immediately.
• Do not touch the lens.
• Do not disassemble the thermometer.
• Basic safety precautions should always be observed, especially when the thermometer is used on or near children and disabled persons.
• This thermometer is not intended to substitute for a consultation with your physician.
• This thermometer and the subject must remain in a stable environment for at least 30 minutes before measuring the temperature.
• When the measured temperature falls within the fever temperature range of >=37.5°C (99.5°F) and <43°C (109.4°F), as indicated by the red LED on display, please consult with your physician immediately.

**Restrictions of use:**

This thermometer is clinically proven to produce accurate temperature measurements. However, please be advised that the accuracy cannot be ensured when the thermometer is not clean. Check that the probe is clean before taking a measurement.

**Intended use:**

Universal Infrared Forehead Thermometer is intended for the intermittent measurement and monitoring of human body temperature from forehead measurement at home. A control measurement using a conventional thermometer is recommended in the following cases:

- 1. If the reading is surprisingly low,
2. For new-born infants up to 100 days old,
3. For children under three years of age who have a weakened immune system or who react unusually in the presence or absence of fever.

**How does it work:**

The thermometer measures the infrared heat generated by the surface of the skin over the vessel and its surrounding tissue. The thermometer then converts it into a temperature value.

Note: The thermometer does not emit any infrared signal.

**DETAILED INFORMATION**

**About Normal Body Temperature & Fever.**

- The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally.
• Vasoconstriction, an effect which constricts the blood vessels and cools the skin, can occur during the early stages of a fever. In this case, the temperature measured by the Universal Infrared Forehead thermometer may be unusually low. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes.
• As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer.
• Body temperature can vary from one individual/person to next. It also varies by location on the body and time of day.
• Please keep in mind that temperatures measured from different sites, even at the same time, should not be directly compared.
• Fever indicates that the body temperature is higher than normal. This symptom may be caused by infection, overdressing or immunization.
• Some people may not experience fever even when they are ill. These include, but are not limited to, infants younger than 3 months old, persons with compromised immune systems, persons taking antibiotics, steroids or antipyretics (aspirin, ibuprofen, acetaminophen), or persons with certain chronic illnesses.
• Please consult your physician when you feel ill even if you do not have fever.

**MAINTENANCE**

**Care & Cleaning**

- The probe is not waterproof. Please wipe with a clean and dry cotton swab to clean the probe on the inside.
• The body of the thermometer is not water-resistant. Never put the thermometer under a running tap or submerge it into water. Use a soft and dry cloth to clean it. Do not use abrasive cleaners.
• Store the thermometer in a cool and dry location. Free from dust and away from direct sunlight.

**Environment: (discarding this product)**

When discarding this product, please return it to your retailer; they'll ensure environment-friendly processing. Exhausted batteries should be brought to your local collection point for small chemical waste.

**DECLARATION OF CONFORMITY**

The declaration of conformity is available at: <http://DOC.hesdo.com/BC-37-DOC.pdf>

**WARRANTY**

For the product you have a warranty of 24 months from the date of purchase. During this period, we guarantee the free repair of defects caused by material and workmanship errors. All this subjected to the final assessment of the importer.

The warranty becomes null and void: In case of improper use, incorrect connections, leaking and/or incorrectly installed batteries (if applicable), use of unoriginal parts or accessories, negligence and in case of defects caused by moisture, fire, flooding, lightning and natural disasters. In case of unauthorised modifications and/or repairs performed by third parties. In case of incorrect transportation of the device without adequate packaging and when the device is not accompanied by this warranty card and proof of purchase. All further liability, especially concerning any consequential damages, is excluded.

Service Help icons, WWW.ALECTO.NL SERVICE@ALECTO.NL, Alecto logo, Hesdo, Australiëlaan 1, 5232 BB, 's-Hertogenbosch, The Netherlands