Liebe Kundin, lieber Kunde!

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




Product number: 377 413

Dear Customer

Mild pain and tension in the back, arms and legs are unpleasant, but can often be improved with massages.

Your new massage band helps you achieve this, especially when it comes to areas of the body that are otherwise difficult to reach.

We hope you enjoy using your new massage band!

Your Tchibo Team

For your safety

Read the safety warnings carefully and only use the product as described in these instructions to avoid accidental injury or damage.

Keep these instructions for future reference.

If you pass this product on to another person, remember to give them these instructions.

Intended use

The massage band is intended for self-massage on the back and legs to increase your sense of well-being in everyday life.

The massage band is intended for household use.

It is not suitable for commercial use in fitness centres or therapeutic institutions.

The massage band is not a toy.

If you allow children to use the massage band, instruct them on how to use it correctly and supervise them while doing so.

Consult your doctor!

- **Before starting a self-massage routine, consult your doctor and ask how much massaging you should be undertaking.** Improper use means that it is not possible to rule out damage to your health.
- If you have specific physical impairments, for example if you wear a pacemaker, have cardiac problems, if you are pregnant, or if you suffer from inflammation of the joints or tendons, water retention

in the arms or legs, orthopaedic complaints, or if you have had surgery, you must arrange an exercise programme in accordance with your doctor.

- Stop massaging immediately and see your doctor if any of the following symptoms occur: nausea, dizziness, excessive shortness of breath or pain in the chest region.
You should also stop the massage immediately if you experience headaches or pain in your joints and muscles.
- If you are susceptible to dizziness, migraines or headaches, you should ask your doctor whether neck massages are good for you.
- In the case of conditions affecting the veins and connective tissues, varicose and spider veins, the leg massages should only be performed using gentle pressure. Discuss the type and extent of the massages with your doctor.
- Not suitable for therapeutic massages!
Not intended for lymphatic massages.

Do not use the massage band...

... if you have a herniated disc in the neck area.

... if you have suffered pathological changes in your vertebrae.

... if you have consumed any medication or alcohol, as this can reduce your perceptive abilities.

... on your throat or in the area of genitals or kidneys.

... if your skin shows signs of swelling, inflammation, injuries, lacerations, broken capillaries, varicose veins, acne, couperose, herpes or other skin infections.

... on babies and small children! Children, disabled people, the frail and people in need of care should only use the massage roller under adult supervision.

Important notes

- Wear comfortable clothing.
- Tie up long hair to ensure that it does not become stuck in the rolls. Do not wear any items of jewellery or clothing with straps, fringes, etc., as these may become stuck in the rolls.
- Always inspect the massage band before use. Do not use it if it shows any signs of damage.
- Do not make any modifications to the product.
- Only use the product for its intended purpose. Do not hang from it, use it to hang any objects or to transport any loads, etc.

What to consider when massaging

- Back massages: Stay upright and balanced - never arch your back!
- Tense your bottom and stomach, but be sure never to make them cramp! The massage is supposed to relax you. Listen to your body! Move the massage band over your muscles in such a way that the pressure feels good.
- If your neck area is particularly sensitive, or if your neck muscles are weak, massage yourself gently and only for as long as it feels pleasant and soothing.
- A massage should last no longer than 10 minutes on any single area of the body. Massaging your neck too firmly or for too long can lead to tension instead of the desired relaxation.
- Move the band over the muscles at a calm, even pace. Keep breathing evenly while doing so.

How to use

- ▷ Hold the massage band firmly by its loops and move it over each body part in gentle motions, but with adequate pressure.

Care

- ▷ If necessary, wipe the massage band down with a damp cloth. Do not use any harsh products, or any sharp or abrasive objects for cleaning.
- ▷ After cleaning, allow the massage band to air-dry at room temperature. Do not place it on a heater and do not use a hairdryer or any similar appliances!
- ▷ Store the massage band in a cool, dry place. Protect it from sunlight as well as any pointed, sharp or rough objects and surfaces.
- ▷ Do not leave the product on sensitive furniture or floor surfaces. Some varnishes, synthetic substances and furniture/floor care products may react adversely with the material of the product and soften it. To avoid any unwanted traces, store the product in a suitable box.

- ▷ Due to manufacturing processes, there may be a slight odour when you first remove the product from its packaging. However, this is completely harmless. Unpack the product and leave it to air out well. The odour will disappear after a short while. Make sure there is sufficient ventilation.

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Tchibo GmbH, Überseering 18, 22297 Hamburg, Germany