
Ultrakompakte Badwaage

de Bedungsanleitung

Tchibo GmbH D-22290 Hamburg · I38557FV03X00XIV · 2023-08

Artikelnummer: 680 609

Liebe Kundin, lieber Kunde!

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
Tchibo GmbH D-22290 Hamburg · I38557FV03X00XIV · 2023-08

Product number: 680 609

 www.tchibo.de/instructions
(Please enter the product number in the box labelled “Bedienungsanleitungssuche” and click on “Suchen”)

Safety warnings

Read the safety warnings carefully and only use the product as described in these instructions to avoid accidental injury or damage. Keep these instructions for future reference. If you give this product to another person, remember to also include these instructions.

Intended use

The scales are designed to weigh people of up to a maximum weight of 150 kg / 330 lb. They are designed for home use and are not suitable for commercial purposes.

Danger to children

- Keep the packaging material out of the reach of children. It poses several risks, including the risk of suffocation!
- Swallowing batteries can be fatal. Swallowing a battery can cause severe internal chemical burns and even death within 2 hours. Therefore, keep both new and old batteries as well as the product itself out of the reach of children. If you think that a battery has been swallowed or has entered the body in any other way, seek medical advice immediately.

Danger of fire/explosion

- Standard batteries must not be charged, dismantled, thrown into fire or short-circuited.

Risk of injury

- Do not stand on the scales if your feet are wet and do not use them if the weighing surface is damp. Danger of slipping!
- The weighing surface is made of glass. The scales must not be dropped and you must not let hard or heavy objects fall onto them. There is a danger that the glass could shatter and cause injury.
- Do not stand on the edge of the scales; they may tip over.
- If a battery has leaked, avoid all contact with skin, eyes and mucous membranes. If necessary, rinse the affected areas with water and consult a doctor immediately.

Material damage


- The scales must not be exposed to excessive moisture. Never immerse them in water!

- Do not store the scales in a hot location. Do not expose them to direct sunlight.
- Do not place any items on the scales except for the purpose of weighing them.
- Do not overload the scales as this could damage the weighing cells. The maximum total weight of 150 kg must not be exceeded!
- Do not subject the product to extreme fluctuations in temperature, moisture, direct sunlight, dust, strong impacts or jolts. This could damage the sensitive measuring technology used in the scales.
- Do not make any modifications to the scales and do not open them. Only have repairs to the scales carried out by a specialist workshop.
- Protect the batteries from excessive heat. Remove the batteries from the product if they are flat or if you will not be using the product for a prolonged period. This helps to avoid damage due to leaking.
- Always replace all batteries at the same time. Do not mix old and new batteries, different types or brands of battery, or batteries with different capacities.
- When inserting the batteries, ensure correct polarity (+/-).
- Clean the contacts on the batteries and the device as required before inserting batteries. Risk of overheating!
- Do not use any caustic, harsh or abrasive products for cleaning.
- Some varnishes, synthetic substances and floor care products may react adversely with the material of the non-slip feet and soften them. If necessary, place a non-slip mat under the product.


Malfunctions

- Keep the scales away from devices that emit strong electromagnetic fields, such as radio-controlled appliances, mobile phones, etc., as they can cause the scales to malfunction, causing incorrect values to be displayed. If the scales have malfunctioned and/or display incorrect values, wait a moment until the scales switch off automatically. Then place the scales in a different location and try using them again.

At a glance (features)


Safety warnings on the protective film

Inserting / replacing batteries

1. Place the scales with the base facing upwards on a soft cloth or something similar to prevent any scratches to the surface.
2. Slide the battery compartment cover in the direction of the arrow to remove it.

3. Insert the batteries as illustrated, starting with the two outer batteries before inserting the middle one. Ensure that the polarity (+/-) is correct.
4. Slide the battery compartment cover back on.

- The first measurement taken after the batteries have been inserted or replaced is used to calibrate the scales - disregard this measurement. **C** appears briefly on the display and then **O.O**. The scales are then ready for use.
- The batteries need to be replaced if **Louu** appears on the display.

Cleaning and storage

- Wipe the scales as needed with a damp cloth. Let the scales dry thoroughly before stepping onto them again. Otherwise you may slip!

Technical specifications

Model:	680 609
Max. weighable weight:	150 kg / 330 lb
Display accuracy:	0.1 kg / 0.2 lb
Units of measurement:	kg / lb
Batteries:	3x LR03 (AAA)/1.5 V (ZnMnO2)
Ambient temperature:	+10 to +40°C
Made exclusively for:	Tchibo GmbH, Überseering 18, 22297 Hamburg, Germany

In the course of product improvement, we reserve the right to make technical and optical modifications to the product.

Use

Setting the unit of measurement

You can choose between two weight units:


- kilograms (**kg**)
- pounds (**lb**)

The default unit of measurement is kilograms (**kg**).

1. Place the scales on a solid, level surface.
2. Push down onto the scales firmly for a moment until the display switches on. However, do not strike the surface of the scales and do not press it erratically.
3. Lift up the scales and press the **UNIT** button on the underside until the desired units are shown on the display.

Weighing

- Place the scales on as level a surface as possible before weighing. This ensures accurate measurements.
- Do not place the scales on a soft surface (e.g. carpets, etc.) as this may lead to inaccurate results. The surface must be hard and even.

- Step onto the scales as illustrated. Stand still and do not hold onto anything.

As soon as the display flashes 3 times, the scales have determined the weight and this is then displayed for approx. 10 more seconds, even if you step off the scales.

- Step off the scales. After approx. 10 seconds, the display will switch off.

Auto-calibration

From time to time the scales calibrate themselves. If this is the case, **C** will be shown on the display after the scales are switched on or after weighing.

- Wait until **O.O** is shown again before standing on the scales.

Tips for obtaining precise measurements

- To obtain comparable results, you should always try to weigh yourself under the same conditions. Therefore, always weigh yourself:
 1. at the same time of day,
 2. before eating,
 3. without clothes, if possible,
 4. in the same place,
 5. at the same ambient temperature, if possible.
- To obtain correct weighing results, the surface must be level and as hard as possible. Soft floors such as carpets, etc. can falsify the result.
- If the air humidity is too high, an incorrect result may be displayed. It is therefore best to weigh yourself before showering or bathing.
- The weight of an adult person can vary by about 2 kg over the course of one day.

Capacity overload

NOTICE - material damage
Do not overload the scales, as this could damage the weighing cells. The maximum total weight of 150 kg must not be exceeded!

If **Err** appears in the display, the maximum weight of 150 kg has been exceeded.

- Step off the scales or remove the weight immediately to prevent damage to the scales.

Resetting the scales

If the scales are clearly displaying incorrect results, reset them to the factory settings:

- Remove the batteries as described in “Inserting / replacing batteries” for a few seconds.

The first measurement taken after the scales have reset is used to calibrate the scales - disregard this measurement. **C** appears briefly on the display and then **O.O**. The scales are then ready for use.

Problems / solutions

Nothing is shown on the display.

- Are the batteries flat?

The scales show an unrealistic result.

- Have the scales been dropped or exposed to severe mechanical shocks/knocks? Reset the scales (see “Resetting the scales” in the chapter “Use”).
- Has the correct unit of measurement been set?
- Are the scales on a suitable (solid, level) surface?
- Is there a mobile phone or other radio-controlled device nearby? This can interfere with the results. Reset the scales (see “Resetting the scales” in the chapter “Use”) and, if necessary, place the scales outside the electro-magnetic field.
- Internal calibration after replacing the battery or resetting the scales: weigh yourself again.

Disposal

The product, its packaging and the supplied batteries have been manufactured from valuable materials that can be recycled. Recycling reduces the amount of refuse and helps to preserve the environment.

Dispose of the **packaging** at a recycling point that sorts materials by type. Make use of the local facilities provided for collecting paper, cardboard and lightweight packaging by type.

 **Appliances, single-use batteries and rechargeable batteries** marked with this symbol must not be disposed of along with household waste!

You are legally obliged to dispose of old appliances separately from household waste. Electronic appliances contain hazardous substances and, if stored or disposed of improperly, may cause harm to health and the environment. Information about collection points where old appliances can be disposed of free of charge is available from your local authorities. Flat single-use and rechargeable batteries must be handed in at a local authority collection point or returned to a battery retailer.